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English PACIFIER THERMOMETER

Please read the instructions in this guide carefully before use.

DESCRIPTION OF EQUIPMENT

CHARACTERISTICS

1. LCD readout to indicate the measured temperature.
2. Recall function for the last 5 measurements stored in the memory.
3. Audible signal indicating that the maximum temperature has been reached.
4. Self-test and automatic shut-off function.

SPECIFICATIONS

| | |
|---------------------|---|
| Temperature range: | 32 to 44°C |
| Tolerance: | Between 35 and 41°C: +/- 0.1°C. Below 35°C and above 41°C: +/- 0.2°C |
| Graduations: | 0.1°C |
| Working conditions: | 0 to 50°C Humidity: 10 to 83 % |
| Storage conditions: | -20 to 60°C Humidity: 10 to 83 % |
| Battery: | 1.55 V DC Alkaline Manganese battery, type SR41, LR41, UCC392, UCC192 |
| Battery life: | Approx. 100 hours for continuous use or approx. 18 months when used for 10 minutes per day. |
| Sensor type: | Thermistor |
| Readout: | 5 mm Liquid Crystal Display (LCD) |
| Memory: | Indication of last temperature measured on LCD readout |
| Alarm: | For approximately 4 seconds once the maximum temperature reached has stabilised |
| Product dimensions: | 55 x 50 x 50 mm |
| Weight: | Approximately 12 grams with battery |

INSTRUCTIONS FOR USE

1. Hold the On/Off button pressed down until the indication "188.8 C/M" followed by 1 M is displayed on the LCD readout. Release the button. The screen displays L°C/M, then 36.5°C. You can use the thermometer as soon as the °C symbol begins to flash. To recall the last 5 measurements, hold the On/Off button pressed down until "188.8 C/M" appears on the LCD readout, followed by the indication 1M, 2M, ..., 5M obtained after each pressing of the button.
2. If the ambient temperature is greater than 44.0°C, the LCD readout displays H°C. If it is less than 32.0°C, the LCD readout displays Lo°C.
3. Place the pacifier bulb in the child's mouth and allow the child to suck on the pacifier thermometer for at least 3 minutes.
4. Once the temperature has stabilised on the maximum measurement for 3 minutes, an audible signal is emitted for 8 seconds.
5. Remove the pacifier thermometer from the child's mouth and read the temperature indicated on the LCD readout.
6. Stop the thermometer by pressing the On/Off button again. To save the battery, we recommend switching the thermometer off as soon as you have noted the temperature measurement. However, the thermometer will switch itself off automatically if no temperature measurement is made within 10 minutes.
7. If you have to interrupt the measurement for any reason, repeat the procedure described above.

PRECAUTIONS FOR USE

1. Before each use, examine the pacifier bulb for cracks, tears, lacerations or other signs of wear. Do not use the pacifier thermometer if it is worn.
2. Do not expose to direct sunlight, high temperatures, moisture or dirt.
3. Do not attempt to pull on the nipple or twist it to remove it from the device as you may damage the temperature sensor and electronic components and thereby impair the sensitivity of the thermometer.

4. Do not allow a child to play with the thermometer. It is not a toy. It is a sensitive electronic device and may be damaged if it is not correctly used under adult supervision.
5. Use this thermometer for oral temperature measurements only.
6. The batteries are not replaceable. Do not attempt to remove it.

CLEANING

The thermometer can be sterilized by immersing ONLY the silicone rubber nipple in a sterilizing solution such as alcohol. Immerse the nipple in the solution for 1 second ONLY. The protective clip-on plastic cap can be used as a container for this purpose.

Warning: The Pacifier Thermometer should never be completely immersed in sterilizing solution or sterilized in boiling water as this will damage the electronic components.

Do not throw old batteries away. Place them in a container specially provided for worn batteries so that they may be collected and recycled.

WARRANTY

This product is guaranteed for 2 years against material and manufacturing defects. During this period, such defects will be repaired free of charge (proof of purchase must be presented for claims under the guarantee). This guarantee does not cover damage resulting from accidents, incorrect use or negligence. For any claims, please first contact the store in which you purchased the product.

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Terraillon SAS France & Headquarters
Terraillon UK Ltd
Terraillon Asia Pacific Ltd
Terraillon Corp USA